



### 1 - NYLON PLASTIC COMPONENTS

#### For cleaning

Use soft cloth and non-abrasive neutral detergent. Do not use alcohol, solvent-based, acid or abrasive cleaners.

#### For disinfection

We recommend the use of sodium hypochlorite (0.1% - 0.5% sol.) or hydrogen peroxide.\*



### 2 - POLYPROPYLENE PLASTIC COMPONENTS

#### For cleaning

Use soft cloth and non-abrasive neutral detergent. Do not use alcohol, solvent-based, acid or abrasive cleaners

#### For disinfection

Denatured alcohol, sodium hypochlorite (0.1% - 0.5% sol.), ethanol (70% sol.), hydrogen peroxide (0.5% sol.) and benzalkonium chloride can be used.\*



### 3 - POLYURETHANE PLASTIC COMPONENTS

#### For cleaning

Use soft cloth and non-abrasive neutral detergent. Do not use alcohol, solvent-based, acid or abrasive cleaners

#### For disinfection

Denatured alcohol or ethyl alcohol in a maximum concentration of 5% is recommended.\*

#### PRODUCT NAME

**BeO 1D**

#### TYPE

**Adjustable armrest**

#### MATERIALS

**Nylon plastic**

**Polypropylene plastic**

**Polyurethane plastic**

**Before disinfecting the entire surface of any component, it is advisable to proceed with tests on a small part of the product.**

\* Analysis of the effects resulting from the application of chemicals on product surfaces, performed by an accredited laboratory.