



NYLON PLASTIC COMPONENTS

For cleaning

Use soft cloth and non-abrasive neutral detergent. Do not use alcohol, solvent-based, acid or abrasive cleaners.

For disinfection

We recommend the use of **sodium hypochlorite** (0.1% - 0.5% sol.) or **hydrogen peroxide**. *



PAINTED ALUMINIUM COMPONENTS

For cleaning

Use soft cloth and non-abrasive neutral detergent. Do not use alcohol, solvent-based, acid or abrasive cleaners

For disinfection

The use of **sodium hypochlorite** (0.1% - 0.5% sol.) is recommended, while the use of **hydrogen peroxide** is not recommended. **Denatured alcohol** and **ethanol** may generate some slight aesthetic changes. *

Before disinfecting the entire surface of any component, it is advisable to proceed with tests on a small part of the product.

PRODUCT NAME

Plane 100x100

Plane Major

TYPE

Fixed support

MATERIALS

Painted aluminium

Nylon plastic

* Analysis of the effects resulting from the application of chemicals on product surfaces, performed by an accredited laboratory.