



#### POLYPROPYLENE PLASTIC COMPONENTS

##### For cleaning

Use soft cloth and non-abrasive neutral detergent. Do not use alcohol, solvent-based, acid or abrasive cleaners

##### For disinfection

Denatured alcohol, sodium hypochlorite (0.1% - 0.5% sol.), ethanol (70% sol.), hydrogen peroxide (0.5% sol.) and benzalkonium chloride can be used.\*

**Before disinfecting the entire surface of any component, it is advisable to proceed with tests on a small part of the product.**

#### PRODUCT NAME

**Modula 1**

#### TYPE

**Seat adapter**

#### MATERIALS

**Polypropylene plastic**

\* Analysis of the effects resulting from the application of chemicals on product surfaces, performed by an accredited laboratory.